

# LCA

## LED CELL ACTIVATOR

Effective for activating cell  
or growth factor before using



LCA  
CELL ACTIVATOR



N-BIOTEK  
Leading Biotechnology



**N-BIOTEK**  
Leading Biotechnology

# LED CELL ACTIVATOR

## is effective for activating cell or growth factor before using.

### FEATURES

- Selective wavelength according to the purpose of treatment.
- RED : to Help cell's growth, collagen synthesis
- BLUE : to Help cell's self-purification, to restrain lipolysis
- GREEN : to Help cellular immunity increase.
- YELLOW : to Stimulate cell's activation.
- N-BIOTEK'S Patent product.

### SPECIFICATION

LCA (LED Cell Activator)	
LED Light	Red/Yellow/Green/Blue
Time	
range	Continuous or up to 47h 59min
accuracy	± 1%
increment	1 Minute
Control	Digital PID controller
Display	5 Digit LEDx2ea
Operating Panel	Touch Button
Block Material	Solid anodized aluminum
Block Capacity	50ml Conical Tube 1Holex2ea
Dimension (in)	30Øx100mmx2ea
(out)	250(W)x250(D)x175(H)mm
Power	110/220V, 50/60Hz, 30W
Weight	5kg



# N-BIOTEK

Leading Biotechnology

402-803 Techno-Park Yakdae-dong,  
 Wonmi-gu, Bucheon-si, Gyeonggi-do,  
 420-734, KOREA  
 TEL: 82-32-321-2100  
 FAX: 82-32-328-2372  
 E-mail: webmaster@n-biotek.com  
 Web: www.n-biotek.com

CO2 INCUBATORS  
 IR CONCENTRATORS  
 LIVE CELL STATIONS  
 SHAKERS & WATER BATHS  
 WORK STATIONS  
 OTHERS  
 BIO MEDICAL EQUIPMENT  
 BIOLOGICALLY CLEAN ROOM

**N-BIOTEK is a global-top  
 biotechnology leading  
 company to develop,  
 supply, manufacture  
 Beauty, Biomedical,  
 life science equipments**