




RG-03/06/12

#### RG-03/06/12 – RM/RTM, Powerful High Viscosity Mixers

##### Features:

- Stirring capacity maximum 80Liter (H<sub>2</sub>O).
- Suitable for material with high viscosity or large volume.
- Chemical resistant coating as well as stainless steel chuck and shaft.
- Small motor with high torque, quiet and vibration less during working.
- Stirring shaft can be adjusted up & down freely.
- Auto function of current accelerating and restraining.
- Overheat and overload safety device.



##### Options:

- RM type with needle tachometer.
- RTM type with needle tachometer & torque meter.

Model	RG-12	RG-06	RG-03
Motor	DC-200W		
Speed (rpm)	200 ~ 1200	100 ~ 600	50 ~ 300
Torque (kg-cm)	15	30	60
Body size (mm)	W125xD250xH250		
Net weight (kg)	6.5		
Power	AC 110/220V, 50/60Hz		



RC-350RM

#### RC-Series, High Torque DC-Mixers

##### Features:

- Auto function of current accelerating & restraining.
- Small motor with high torque provides quiet and vibration less operation.
- Over heat protection.

##### Option:

- "D" model with digital rpm meter.
- RM type with needle tachometer & Ammeter.



##### Accessories:

RC-450RM/D Shaft: Ø8x500mm

RC-200RM/D RC-350RM/D Shaft:

Max. Ø13x280mm

Propeller: No. 11 Ø50mm, Ø80mm.

Stand base: W400xD315xH54mm, 8kg.

Holder

Pole

494

Model	RC-200 RM/D	RC-350 RM/D	RC-450 RM/D
Motor	DC 200W	DC 350W	DC 450W
Speed (rpm)	4,000	7,000	4,000
Torque (kg-cm)	7	5.4	18
Net weight (kg)	17	18	25
Power	AC 110/220V, 50/60Hz		AC 220V, 50/60Hz