

## Multichannel Pipette

- Manifold can rotate 360° to allow any convenient pipetting angle
- Specially designed tip cone with O-ring resulting in:
  - Leak free operation
  - Smooth as well as light loading and ejection
  - Choice of variety of tips
- Each channel is calibrated individually
- Includes all innovative features of Prime single channel pipette

### 8 Channel Micropipette

Cat. No.	Volume Range [ µl ]	Increments [ µl ]
SMA8-010	0.5 - 10 µl	0.02
SMA8-020	2 - 20 µl	0.02
SMA8-050	5 - 50 µl	0.1
SMA8-100	10 - 100 µl	0.1
SMA8-200	20 - 200 µl	0.2
SMA8-300	30 - 300 µl	0.2



### 12 Channel Micropipette

Cat. No.	Volume Range [ µl ]	Increments [ µl ]
SMA12-010	0.5 - 10 µl	0.02
SMA12-020	2 - 20 µl	0.02
SMA12-050	5 - 50 µl	0.1
SMA12-100	10 - 100 µl	0.1
SMA12-200	20 - 200 µl	0.2
SMA12-300	30 - 300 µl	0.2

### Variable Volume Micropipette

Cat. No.	Volume Range [ µl ]	Increments [ µl ]
SVA-125	0.1-2.5	0.002
SVA-100	0.5-10	0.02
SVA-200	5-50	0.1
SVA-500	10-100	0.1
SVA-600	100-1000	1
SVA-700	2-20	0.02
SVA-800	20-200	0.2
SVA-905	500-5000	10
SVA-1000	1000-10000	20

### Fix Volume Micropipette


Cat. No.	Volume [ µl ]
SFA-205	2.5
SFA-5	5
SFA-10	10
SFA-20	20
SFA-25	25
SFA-50	50
SFA-100	100
SFA-200	200
SFA-250	250
SFA-500	500
SFA-1000	1000
SFA-2000	2000
SFA-5000	5000
SFA-10000	10000

### Pipette Stand

Cat. No.	Description
SAS-01	Carousel stand for 6 pipettes
SAS-02	Shelf mount pipette holder



# Where Form Meets Function



## Fine Care Corporation

194/6/2, Dantali Industrial Estate, Gota-Vadsar Road,  
Dantali, Dist.: Gandhinagar - 382 721, Gujarat, India.  
Phone +91-2764-286641 Fax : +91-2764-286643  
E-mail : info@pfact.in Web : www.pfact.in

