

DG-600, Personal Ultrasonic Cleaner (without heating)



| Model | DG-600 |
|-------------------------|---------------------------|
| Operating frequency | 43 KHz |
| Tank Capacity | 600 ml |
| Internal tank size (mm) | L170xW90xH55 |
| External dimen. (mm) | L193xW113xH110 |
| Internal tank stuff | SUS304 # 0.7t |
| External stuff | ABS with fireproof |
| Dust cover stuff | Semi transparent plastic |
| Timer | 5 min |
| Ultrasonic power | 50W |
| Heateing Power | Without |
| Net weight | 1 kg |
| Power source | AC100~120V/AC220V 50/60Hz |

DG-1300, Personal Ultrasonic Cleaner (with heating)



| Model | DG-1300 |
|-------------------------|--|
| Operating frequency | 42 KHz |
| Tank Capacity | 1,300 ml |
| Internal tank size (mm) | L175xW145xH70 |
| External dimen. (mm) | L228xW183xH158 |
| Internal tank stuff | Stainless steel SUS304 |
| External stuff | ABS with fireproof |
| Dust cover stuff | Semi transparent plastic |
| Time Control | 5 Cycles with Digital time display |
| Ultrasonic power | 60W |
| Heateing Power | 60°C |
| Net weight | 1.75 kg |
| Power source | AC 110-120V 60Hz AC 220-240V 50/60 Hz |

Features:

- High-power ultrasonic transducer for superior cleaning effect.
- Large capacity tank for professional cleaning.
- Add heating for special cleaning demands.
- 3-59 minutes adjustable with digital display.
- 5 Cycles adjustable with digital display.
- Memory function for the power failure.
- Transparent lid with blue light for view the cleaning process.
- With radiator for longer working time.
- Secure circuitry design with moisture-proof PCBA.

Electron workshop, Glasses shop, Jewelry shop, Mobile phone, maintain shop, Hospital Dentistry, Office and household articles etc.

CLEAN-01, 4.5liter Ultrasonic Cleaner (Dental)



- Cycle time set: 2-30 minutes.
- Temp set: 20-80°C.
- Safety feature: Heating stops when tank is empty.
- Water outlet is installed to drain liquid from the tank easily.
- Low noise fits clinical environment.
- Cleaning basket holds delicate instruments in place.
- Keypads & digital display are user friendly.
- Rinse time can be set from 2 to 30 min, rinse temperature can also be set. Power is automatically cut off when is lack of water.
- External drainage system makes easy change of rinse liquid; Low noise design makes quiet operation.
- Specialized basket is designed for delicate minute devices that should be free from friction and braise, for example, high-speed and low-speed hand piece turbines
- Clear operation plate with digital display makes it easy to operate.