



PRINTER-1RG



TM-U220A



CLP521

PRINTER Series

Printer "Kafka" enables two-way thermal printing by mobile head. The design of the mark: matrix 8 × 8 points. Print Speed: 0.75 line / sec.

PRINTER 1/RG

base functions, (main current supply).

PRINTER 1/Z/RG

real time clock, (main current supply).

PRINTER 1/RG Portable

base functions, (battery power and main current supply).

PRINTER 1/Z/RG Portable

real time clock, (battery power and main current supply).

PRINTER SQS

date, time and statistic functions, (main current supply).

TM-U220 Series

TM-U220A

with paper roll reeler & paper cutter.

TM-U220B

with paper cutter.

TM-U220D

basic version.

CL Series

CLP521

Print: thermal, Head: 200dpi, Print Speed: 100mm/sec.

CLP621

Print: thermal, thermo-transfer, Head: 200dpi, Print Speed: 100mm/sec..

CLP631

Print: thermal, thermo-transfer, Head: 300dpi, Print Speed: 175mm/sec..

CL-S700

Print: thermal, thermo-transfer, Head: 200dpi, Print Speed: 250mm/sec.

Model	Microbalances	Analytical balances	Precision balances	Moisture analyzers
Anti-vibration table - mild steel construction	•	•	•	•
Anti-vibration table - stainless steel construction	•	•	•	•
Laboratory weighing table	•	•		
Additional LCD display WD-3		•	•	
Anti-draught shield			•	
Anti-static ionizer	•	•	•	
Density determination kit		•	•	
Environmental module	•	•	•	
Foot button for Tare and Print functions		•	•	
Pipettes calibration kit	•			
Power adapter with battery and charger		•	•	
Rack for under hook weighing		•	•	
Set for the determination of air	•	•		

Model	Microbalances	Analytical balances	Precision balances	Moisture analyzers
Set for determination of water vapour permeability				•
Suitcase for balance			•	
Computer software Pipettes	•	•		
Computer software	•	•	•	•
Thermal printer 1/RG	•	•	•	•
Thermal printer 1/RG Portable	•	•	•	•
Thermal printer 1/Z/RG		•	•	•
Thermal printer 1/Z/RG Portable		•	•	•
Thermal printer SQS		•	•	
Needle printer TM-U220A	•	•	•	
Needle printer TM-U220B	•	•	•	
Needle printer TM-U220D	•	•	•	
Label printer CLP 521		•	•	
Label printer CLP 621		•	•	
Label printer CLP 631		•	•	
Label printer CLP S700		•	•	
Disposable pans (set of 50 pcs.)				•
Infra-red emitter				•